



SkinPen<sup>®</sup>  
BY CROWN AESTHETICS

Uncover your best skin

SKINPEN<sup>®</sup> PRECISION.




**Advanced Cartridge Unit (ACU) contains reciprocating mechanism that ensures long-term reliability**

14 medical-grade, stainless steel needles validated to remain sharp for at least 80 minutes of use on scar tissue

100% ultrasonic sealed cartridge

Proven precision depth dial penetrates at set depth

Single-use lockout feature prevents accidental reuse



# World-Class Design and Technology

LED smart technology ensures optimized procedure

Cordless design allows maximum versatility during procedure

Inductive charging ensures full charge for every procedure



## Superior Features. Unsurpassed Safety.

**U.S.A. ENGINEERED AND MANUFACTURED**

SkinPen® Precision is backed by years of multiple peer-reviewed publications and 90 independent validation studies.<sup>1</sup>



- U.S. FDA-cleared, CE-marked and ISO-certified
- Registered with Health Canada, Australian Register of Therapeutic Goods and New Zealand Medicines and Medical Devices Safety Authority
- Available in 35+ countries
- 500,000+ treatments performed
- More than 6,000 devices sold
- 96% of SkinPen Precision providers expressed high satisfaction<sup>2</sup>

# Improve Fine Lines and Wrinkles

SkinPen® Precision improves fine lines and wrinkles while remodeling skin to look younger.



BEFORE



AFTER FOUR TREATMENTS  
*Results may vary*



BEFORE



AFTER FIVE TREATMENTS  
*Results may vary*



BEFORE



AFTER FOUR TREATMENTS  
*Results may vary*

# Reduce Acne Scars

In clinical studies, both patients and providers saw significant improvement in their skin's overall appearance.<sup>1</sup>



BEFORE



AFTER FOUR TREATMENTS  
*Results may vary*



BEFORE



AFTER SIX TREATMENTS  
*Results may vary*



BEFORE



AFTER THREE TREATMENTS  
*Results may vary*



BEFORE



AFTER TWO TREATMENTS  
*Results may vary*

## Post-Procedure Reactions

Some patients have experience post-procedure reactions that include dryness, roughness, tightness, itching, peeling, discomfort and tenderness lasting between one and six days after the procedure.<sup>1,3</sup>

# Patients Love SkinPen® Precision

## % Patients, Six Months Post-Treatment




**85%** of patients stated that they were satisfied with SkinPen Precision

**90%** of clinical trial patients would recommend SkinPen Precision to others.<sup>1</sup>

## Longterm Results

The percentage of SkinPen clinical trial subjects reporting improvement in the appearance of facial acne scars over time.



Results from a clinical study of 20 patients treated with three SkinPen Precision treatments, spaced one month apart. Assessment was made at six months. Provider assessment was made using the Clinician's Global Aesthetic Improvement Scale and patient assessments were made using the Subject Global Aesthetic Improvement Scale.<sup>1</sup>

SkinPen Precision is a precise micromechanical skin remodeling device that initiates skin remodeling to treat scar tissue including acne and surgical scars, and may also be used for general cosmesis rejuvenation to improve fine lines, wrinkles, pigmentation and overall skin appearance.





# SKINPEN GIVES YOUR SKIN THE RESET THAT IT NEEDS IN AS FEW AS THREE TREATMENTS\*

Hitting reset on your skin can be as simple as a 30-minute microneedling treatment. SkinPen creates micro-channels (microscopic holes) at a rate of 1600 per second. These wounds stimulate your skin, kicking its healing process into high gear. Best of all, SkinPen doesn't rely on heat or chemicals, so the results are natural and safe for all skin tones and types.



## POPULAR

- 90% of clinical study patients would recommend SkinPen Precision to family and friends
- Recommendation of three to six treatments ensures patients will return
- Happy patients are more likely to return to your practice



## AFFORDABLE

- Cost-effective treatment to expand your patient base
- High ROI for your practice
- Gateway treatment brings in new aesthetic patients



## MINIMALLY INVASIVE

- Little to no downtime
- Most patients get back to their regular routine right away
- Attractive to patients with busy schedules



## TIME-EFFICIENT

- Short appointments
- Patients can be in and out in 30-45 minutes
- Fit many patients into your schedule



## DELEGABLE

- Easy treatment
- May be performed by trained and certified clinical staff



## SUPPORTED

- Ongoing training support
- Marketing materials available for local distribution



SkinPen® Precision, the superior choice for precise micromechanical skin remodeling.

- Clinically proven to improve acne scars
- Expands your patient base with gateway treatment
- Exceptional value for your practice and patients

References: 1. Data on file. Bellus Medical. 2017. 2. Data on file. Bellus Medical. 2019. 3. SkinPen Instructions for Use.

SkinPen Precision is a precise micromechanical skin remodeling device that initiates skin remodeling to treat scar tissue including acne and surgical scars, and may also be used for general cosmesis rejuvenation to improve fine lines, wrinkles, pigmentation and overall skin appearance. Manufactured by Bellus Medical. For international distribution only.